



## Restriction Enzyme Digest & Gel Electrophoresis Inquiry

Partner Names (& Initials): \_\_\_\_\_  
\_\_\_\_\_

For this experiment, you and your laboratory group will perform two procedures as outlined below:

- The **control procedure** is a method for digesting lambda bacteriophage DNA with restriction enzymes.

As your team performs the control procedure:

- Read each step carefully.
  - Pay attention to how each step is done and think about how that particular step effects all of the other steps, and the outcome of the procedure.
  - Share your observations and questions with your lab team. This should help your team to identify something to change in the experimental procedure that you do next.
- The **experimental procedure** is a method that you will perform that is exactly like the control procedure except that you will **change** one detail. You can choose to change any single detail that appeals to you. You could alter the amount or alter the type of an ingredient in the digestion, such as:
    - Use different restriction enzymes
    - Use a different digestion temperature
    - Use a different length of digestion time

## Restriction Enzyme Digest & Gel Electrophoresis Inquiry

### Control Procedure

Partner Initials: \_\_\_\_\_

### I. Reaction Set-up

Procedure
1. Each group should label 5 clear plastic tubes with their initials and an indicator for the enzyme used: <b>E, C, H or EH</b> ; or <b>(-)</b> for negative control. Put tubes on ice.
2. Use a microcentrifuge to spin down all six colored tubes in your ice bucket. Spin these tubes for 3 to 5 seconds, returning the tubes to the ice bucket when done. <b><u>Be sure that your tubes are BALANCED in the centrifuge!!!</u></b>
3. Add ingredients <b>IN ORDER</b> , according to the following table.

Check off ingredients as you add them:

Reagent	Tube E	Tube H	Tube EH	Tube C	Tube (-)
<b>Nuclease Free H<sub>2</sub>O</b>	3 $\mu$ l	3 $\mu$ l	3 $\mu$ l	3 $\mu$ l	4 $\mu$ l
<b>10x reaction buffer B</b>	1 $\mu$ l	1 $\mu$ l	1 $\mu$ l	1 $\mu$ l	1 $\mu$ l
<b>Lambda DNA (60ng/<math>\mu</math>l)</b>	5 $\mu$ l	5 $\mu$ l	5 $\mu$ l	5 $\mu$ l	5 $\mu$ l
<b>Enzyme</b>	1 $\mu$ l <b>E</b>	1 $\mu$ l <b>H</b>	1 $\mu$ l <b>EH</b>	1 $\mu$ l <b>C</b>	0 $\mu$ l
<b>Total volume</b>	10 $\mu$ l	10 $\mu$ l	10 $\mu$ l	10 $\mu$ l	10 $\mu$ l

4. Spin the reaction tubes for 5 seconds in a microcentrifuge tube to bring the liquid down to the bottom of the tube.
5. Incubate the reactions and the negative control at 37°C for 15 minutes. Record the start time and the total incubation time below: <i>Start Time:</i> _____ <i>End Time:</i> _____ <i>Total Incubation Time:</i> _____ minutes

## Restriction Enzyme Digest & Gel Electrophoresis Inquiry

### *Control Procedure*

Partner Initials: \_\_\_\_\_

### II. Prepare Gel Box with Buffer and Gel

Procedure
1. Obtain a 0.8% agarose gel from your instructor.
2. Loosen the screws and lower the dams on both sides of the gel casting tray.
3. Place the agarose gel tray into the electrophoresis cell with the comb closest to the negative (black) electrode. It should slide in easily. Ensure that the gel tray dams are still down.
4. Add enough buffer to the electrophoresis chamber to just cover the gel, about 325-350ml.
5. Rock the comb very gently, back and forth, in the gel. Gently remove the comb. Set the comb on top of the power supply (please <b>do not</b> throw the comb away).

### III. Load Gel and Run

1. Retrieve your samples from the heat block.
2. Add 2µl of the EZ-Vision DNA stain to each reaction and the negative control.
3. Fill in the blanks (-----) in the loading diagram below, listing the order of the samples loaded on the gel.

With the wells at the top of the gel, lane 1 is the well on the left side.

Lane 1	Lane 2	Lane 3	Lane 4	Lane 5	Lane 6
10µl <b>Negative Control</b> Uncut Lambda DNA	10µl	10µl	10µl	10µl EcoRI/ HindIII Lambda <b>Marker</b>	10 µl
	-----	-----	-----		-----

4. Load 10ul <b>negative control</b> , 10ul prepared EcoRI/HindIII lambda <b>marker</b> ), and load the 10ul of each digest into the appropriate wells as described in the table above.
---

## Restriction Enzyme Digest & Gel Electrophoresis Inquiry

### *Control Procedure*

Partner Initials: \_\_\_\_\_

### **III. Load Gel and Run (continued)**

<b>Procedure</b>
5. Place the lid on the gel box and connect the electrodes into the power supply. Make sure that the black plug is in the black outlet and the red into the red.
6. Turn on the power supply and set it at 250V. Bubbles at the electrodes indicate that electric current is running through the gel. The gel should run for approximately 20 minutes. Record the start time and the total run time below:  <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <span><i>Start Time:</i> _____</span> <span><i>End Time:</i> _____</span> </div> <div style="text-align: center; margin-top: 10px;"> <span><i>Total Run Time:</i> _____ <i>minutes</i></span> </div>

### **IV. Results & Observations**

1. Tape a picture of each of your gels or draw your gel results below. Be sure to annotate your gel picture or diagram with the following information: <ul style="list-style-type: none"> <li>✓ Agarose gel %</li> <li>✓ Name of the running buffer</li> <li>✓ A legend detailing the wells/lanes and you put in them</li> <li>✓ Type of DNA stain (EZ-Vision)</li> <li>✓ Run Time in minutes</li> <li>✓ Voltage at which the gel was run</li> </ul>
2. Record any observations from your results that may be helpful when you go on to perform the experimental procedure.

## Restriction Enzyme Digest & Gel Electrophoresis Inquiry

### *Experimental Procedure*

Partner Initials: \_\_\_\_\_

After performing the control procedure, reflect on how what you did affected your results. Talk to your laboratory team to identify a detail of the control procedure that you will change in the experimental procedure. An **experimental procedure** is exactly like the control procedure except that you will **change** one detail.

Design an experiment according to the following rules:

- Change **only one variable** between your control and experimental procedures. You can choose to change any single detail that appeals to you. Your group may decide to alter the amount or alter the type of an ingredient in the digestion, such as:
  - Use a different restriction enzymes
  - Use a different digestion temperature
  - Use a different length of digestion time
  
- You will be limited to **seven** samples including your control samples.

1. Write down the experimental variable that you've chosen: \_\_\_\_\_

2. Describe the negative and positive controls that you will prepare given the variable you've chosen:

Positive Control(s)	Negative Control(s)

As you and your partner proceed to run the experiment, think about the following:

- **Follow the instructions** to the best of your ability.
- **Write down** any observations that you or your partner make, such as anything that is different than, or in addition to, the instructions as they are written.

# Restriction Enzyme Digest & Gel Electrophoresis Inquiry

## *Experimental Procedure*

Partner Initials: \_\_\_\_\_

### I. Reaction Set-up

Procedure	Observations
1. Write down the experimental variable you have chosen.	
2. Describe the controls that you will run for the experimental variable you have chosen.  <u>Positive (+):</u>  <u>Negative (-):</u>	
3. Devise a way to label your tubes so that initials or symbols can represent what will be in each tube. Record your labelling scheme below:	

<i>Tube Label</i>							
<i>Tube contents</i>							

4. Use a microcentrifuge to spin down all the tubes you will use from your ice bucket. Spin tubes for 5 seconds, returning them to the ice bucket when done. **Remember to BALANCE your centrifuge!!!**

# Restriction Enzyme Digest & Gel Electrophoresis Inquiry

## Experimental Procedure

Partner Initials: \_\_\_\_\_

### I. Reaction Set-up (continued)

Procedure	Observations
5. Below is a recipe for a single 10 $\mu$ l restriction enzyme digest. If you use the recipe for some or all of your reactions, then circle the recipe and annotate it appropriately. If you modify the recipe, record in the open space what your adjusted recipe(s) are.	

<p> <b>3<math>\mu</math>l Nuclease Free H<sub>2</sub>O</b>  <b>1<math>\mu</math>l 10x reaction buffer B</b>  <b>5<math>\mu</math>l Lambda DNA (60ng/<math>\mu</math>l)</b>  <b>1<math>\mu</math>l Enzyme(s)</b> (or water for (-) control)            -----  <b>10<math>\mu</math>l total volume (300ng DNA)</b> </p>	
---	--

6. Spin the reaction tubes for 5 seconds in a microcentrifuge tube to bring the liquid down to the bottom of the tube.
7. Incubate the reactions and the controls. Record the temperature and time of incubation below:  <p style="text-align: center;"> <i>Temperature:</i> _____ °C   <i>Start Time:</i> _____                      <i>End Time:</i> _____   <i>Total Incubation Time:</i> _____ minutes         </p>

## Restriction Enzyme Digest & Gel Electrophoresis Inquiry

### *Experimental Procedure*

Partner Initials: \_\_\_\_\_

### II. Prepare Gel Box with Buffer and Gel

Procedure
1. Obtain a 0.8% agarose gel from your instructor.
2. Loosen the screws and lower the dams on both sides of the gel casting tray.
3. Place the agarose gel tray into the electrophoresis cell with the comb closest to the negative (black) electrode. It should slide in easily. Ensure that the gel tray dams are still down.
4. Add enough buffer to the electrophoresis chamber to just cover the gel, about 325-350ml.
5. Rock the comb very gently, back and forth, in the gel. Gently remove the comb. Set the comb on top of the power supply (please <b>do not</b> throw the comb away).

### III. Load Gel and Run

Procedure	Observations
1. Retrieve your samples from the heat block.	
2. Add EZ-Vision DNA stain to each reaction and the controls at a proportion of 2µl EZ-Vision / 10µl of reaction.	
3. Fill in the blanks (-----) in the loading diagram below, listing the order of the samples loaded on the gel.	

Lane 1	Lane 2	Lane 3	Lane 4	Lane 5	Lane 6	Lane 7	Lane 8
10µl	10µl	10µl	10µl	10µl	10 µl	10 µl	10 µl
-----	-----	-----	-----	-----	-----	-----	-----

4. Remember to reserve space in your gel for one well to hold 10µl EcoRI/HindIII Lambda DNA Marker for a size standard, and load 10µl of your remaining samples into the gel as described in the table above.	
---	--

## Restriction Enzyme Digest & Gel Electrophoresis Inquiry

*Experimental Procedure*

Partner Initials: \_\_\_\_\_

### **III. Load Gel and Run (continued)**

<b>Procedure</b>
5. Place the lid on the gel box and connect the electrodes into the power supply. Make sure that the black plug is in the black outlet and the red in the red.
6. Turn on the power supply and set it at 250V. Bubbles at the electrodes indicate that electric current is running through the gel. The gel will run for approximately 20 minutes.
<p><i>Start Time:</i> _____ <i>End Time:</i> _____</p> <p><i>Total Run Time:</i> _____ <i>minutes</i></p>

### **IV. Results & Observations**

Tape a picture of each of your gel(s) or draw your gel results below or on a separate page and annotate appropriately.

# Restriction Enzyme Digest & Gel Electrophoresis Inquiry

*Experimental Procedure*

Partner Initials: \_\_\_\_\_

## **V. Analysis**

Procedure
1. Describe the differences between how you set up your experimental samples and your control samples.
2. Describe all of the differences that you observed between your control and experimental <u>procedures</u> . This should include the experimental variable or change that you selected, plus any other differences that you observed.
3. Did the gel look the way you expected? Explain why or why not.



**BioPharmaceutical Technology Center**

I N S T I T U T E

# Restriction Enzyme Digest & Gel Electrophoresis Inquiry

*Experimental Procedure*

Partner Initials: \_\_\_\_\_

## **V. Analysis** (continued)

Procedure
4. Based on the <u>results</u> of your experiment, describe how your chosen experimental variable affected your result.
5. Explain what part of your experiment surprised you the most.
6. If you could run a new experiment, what would you try next?